

# dry block heater



## Blocks Dry Block Heater Blocks

A range of standard blocks are available for use with DBH10, DBH20D, DBH30 & DBH30D, including blocks with tapered holes to suit 1.5ml and 0.5ml Eppendorf tubes. Blocks can be manufactured to order for most tubes and small flasks, to suit customer requirements.



[Details / Specifications](#)

### Details

<b>Ratek Product Code</b>	DBH-Blocks
<b>Description</b>	Dry block heater aluminium blocks to suit DBH10, DBH20D, DBH30 & DBH30D
<b>Packaging Contents</b>	Block only

### Features

- Black anodised aluminium
- Temperature stability
- Can be manufactured to order
- Thermometer hole for accurate temperature control

### Warranty period

1 Year

### Download this page in PDF

<ratek-dbhblocks.pdf>

### Specifications

<b>Material</b>	Black anodised aluminium
<b>Heating well layout</b>	Depending on hole size
<b>Overall dimensions</b>	95 x 75 x 50mm
<b>Weight</b>	500 - 1000 grams each depending upon holes

Copyright © 2000 Ratek Instruments Pty. Ltd.

# dry block heater

# dry block heater