

dry block heater



DBH10 Dry Block Heater

A compact block heater to accommodate one standard Ratek test tube block. The solid state proportional temperature control has a range of approx 28°C (ambient + 5°C) to 110°C with a stability of 0.1°C. Each block has provision for a reference thermometer to ensure precise temperature setting.



[Details / Specifications](#)

Details

Ratek Product Code	DBH10
Description	Dry block heater
Packaging Contents	DBH10 dry block heater, instruction sheet

Specifications

Power requirements	240vac/165 watts
Control type	Solid state proportional
Control range	Ambient + 5°C to 110°C (min 28°C)
Safety features	Over temperature cut out
Overall dimensions	W165 x D175 x H90mm
Weight	1.6kg (without block)

Features

- Single-block capacity
- Wide temperature control range
- Over-temperature cutout
- Element indicator

Warranty period

1 Year

Download this page in PDF

[ratek-dbh10.pdf](#)

Copyright © 2000 Ratek Instruments Pty. Ltd.

dry block heater

dry block heater