



BTR10

Blood Tube Roller

Various types of tubes and vials can be accommodated, with a maximum weight of 1Kg per pair of rollers. The BTR10 can hold up to 54 x 5ml blood tubes. The coated aluminium rollers are mounted on eccentric drive spindles to provide a rocking as well as rolling motion, ensuring a thorough mixing action. Parallel rollers can be specified, if preferred. Speed of the rollers is variable between 10 and 60 rpm. To cover a wide range of mixing requirements.

The sprocket and chain drive, coupled with ball race roller mounting, will ensure smooth trouble free operation. The unit can run continuously for many hours unattended.



[Details / Specifications](#)

Details

Ratek Product Code	BTR10
Description	Blood tube roller
Packaging Contents	BTR10 blood tube roller, instruction guide

Specifications

Power requirements	240 VAC/20 watts
Control type	Solid state with variable speed between 10 and 60rpm
Motor	12 volt DC geared
Drive	Chain and sprockets
Rollers	30mm diameter x 300mm long
Overall dimensions	L500 x W400 x H120
Weight	6.5 kg

Features

- Mixes up to 54 x 5ml blood tubes
- Rocking & rolling motion
- Variable speed between 10 and 60 rpm
- Ball race roller mounting for quiet operation

Warranty period

1 Year

Download this page in PDF
[ratek-btr10.pdf](#)