



RPM5

Rocking Platform Mixer

This rocking platform mixer has a variable tilt action to suit a range of mixing and washing applications, with a speed range of 4-57 RPM. The microprocessor controlled drive system is complete with digital tachometer and timer to ensure reliable and repeatable procedures.



[Details / Specifications](#)

General

Included are 3 user preset programs that can be set for speed and run time. The 350mm square platform is mounted on a central pivot, and is fitted with a non slip mat that provides adequate grip for most plates and flat bottom containers. Racks are available to suit most flasks and containers if required.

Details

Ratek Product Code	RPM5
Description	RPM5 Rocking Platform Mixer
Packaging Contents	RPM5 Rocking Platform Mixer, instruction guide

Specifications

Power requirements	240vac/100 watts
Capacity	350 x 350mm platform maximum capacity 6kg's
Mixer control	Solid state with optical feedback
Motor	Brushless ac
Tilt angle	Variable
Speed range	4 - 57 rpm
Timer	0-99 minutes or 0-99 hours with 3 presets
Overall dimensions	W350 x D400 x H200 mm
Weight	10.1kg

Features

- Speed range of 4-57 RPM
- Microprocessor controlled drive system
- 3 user preset programs
- Silent operation
- Tilt-adjustable

Warranty period

1 Year

Download this page in PDF
[ratek-rpm5.pdf](#)