



IT16

Insulated Tank

16-litre insulated water tank for use with the Ratek range of Thermoregulators. Other sizes can be manufactured to customer specifications. The inner tank is fabricated from 304 grade polished stainless steel, and the outer case from a chemical resistant plastic. 25mm of insulation is used between the tank and case to minimise heat loss, allowing for stable temperature control. Gabled or flat lids are also available. It is recommended that tanks to be used with a Thermoregulator are ordered with a bridge plate to protect the Thermoregulator from condensation and minimise heat loss.



[Details / Specifications](#)

Details

Ratek Product Code	IT16
Description	Insulated tank
Packaging Contents	IT16 Insulated tank

Specifications

Tank capacity	Approx 16 litres
Inner tank dimensions	L400 x W200 x D200
External dimensions	L470 x W270 x D245

Features

- Chemical resistant plastic outer
- Even heat transfer
- Suitable for use with TH5 or TH3P thermoregulators
- Bridge plates and gabled or flat lids available
- Drain taps and plumbing fittings can be fitted if required

Warranty period

1 Year

Download this page in PDF
[ratek-it16.pdf](#)

Copyright © 2007 Ratek Instruments Pty. Ltd.