



RT7

Insulated Tank

Insulated water tank for use with the Ratek range of Thermoregulators. Includes plumbing fittings enabling a supply of temperature controlled water to be circulated to an external sealed system. The inner tank is fabricated from 304 grade polished stainless steel, and the outer case from a chemical resistant plastic. 25mm of insulation is used between the tank and case to minimise heat loss, allowing for stable temperature control.



[Details / Specifications](#)

Details

Ratek Product Code	RT7
Description	Recirculating tank
Packaging Contents	RT7 Tank only, user's guide

Specifications

Tank capacity	Approx 7 litres
Inner tank dimensions	L200 x W175 x D200
External dimensions	L270 x W244 x D245
Fittings	2 x Ø12.7mm hose tails
Lid	Stainless steel bidge and flat lid

Features

- 304 BA grade stainless steel used for tank and cover
- Chemical resistant plastic outer
- Even heat transfer
- Suitable for use with TH5 or TH3P thermoregulators

Warranty period

1 Year

Download this page in PDF

[ratek-rt7.pdf](#)

Copyright © 2007 Ratek Instruments Pty. Ltd.